

# M3A Vital Signs Monitor



Outpatient

Transport

Ward



Portable Design



Infrared Ear Thermometer



**EDAN**

A world of potential

# M3A Vital Signs Monitor

EDAN M3A Vital Signs Monitor provides essential measurement of patient's SpO<sub>2</sub>, NIBP and Temperature. It offers an effective solution and brings high performance vital signs monitoring at your finger tips.

<b>30,000</b> NIBP Measurement Review	<b>72 h</b> Trend Review	<b>800</b> Alarm Review
--	-----------------------------	----------------------------

- 3.5" high resolution color TFT-LCD screen & LED display
- Lightweight & portable design
- User-friendly interface
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 17 hours of continuous work
- Compatible with Patient Care Viewer PC data management software
- Nurse call
- Suitable for adult, pediatric and neonatal patients
- Barcode scanner support (optional)
- Ethernet connection for bi-directional communication with MFM-CMS central monitoring system (optional)
- Wi-Fi support with optional AP



## Mounting Solutions



Wall Mount



Rolling Stand



Bed Rail Clamp

## Specification

Configuration	SpO <sub>2</sub>	NIBP	Quick Temp	Infrared Ear Thermometer
M3A(SpO <sub>2</sub> only)	✓	—	—	—
M3A(NIBP only)	—	✓	—	—
M3A(SpO <sub>2</sub> +NIBP)	✓	✓	—	—
M3A(NIBP+Quick Temp)	—	✓	✓	—
M3A(NIBP+Infrared Ear Thermometer)	—	✓	—	✓
M3A(SpO <sub>2</sub> +NIBP+Quick Temp)	✓	✓	✓	—
M3A(SpO <sub>2</sub> +NIBP+Infrared Ear Thermometer)	✓	✓	—	✓



Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,  
1019# Nanhai Rd., Shenzhen | 518067 P.R. China  
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn



© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.

ENG-PM-M3A-V2.0-20150323