



iM8 VET series
Veterinary Monitor

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- In-Vitro Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



Edan Instruments, Inc. | 3/F-B, Nanshan Medical Equipments Park,
1019# Nanhai Rd., Shenzhen | 518067 P.R. China
+86.755.26898326 | www.edan.com.cn | info@edan.com.cn

© 2015 Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission.



ENG-PM-iM8 VET-V2.0-20150323


iM8 VET series

Veterinary Monitor



iM8 VET series

Veterinary Monitor



iM8 VET
12.1" Color TFT-LCD Screen

iM8B VET
10.1" Color TFT-LCD Screen

- High resolution color TFT-LCD Screen
- Defibrillation and electrosurgical interference protection
- Pitch tone (Pulse-tone modulation)
- VGA output(optional)
- Suitable for feline, canine and other animals

Standard Parameters: 3/5-lead ECG, RESP, EDAN SpO₂, NIBP, PR, 2-TEMP

Optional: 2-IBP, Respironics CO₂, EDAN CO₂, Thermal Recorder




Respironics CO₂ (Mainstream/Sidestream, optional)



- Proven ruggedness - Solid state design
- No need to calibrate on a regular basis
- Sidestream sampling rate of 50 ml/min

Display Modes



3/5-lead ECG Analysis

Large Font

Trend Screen

500

NIBP Measurement Review

96 h

Trend Review

60

Alarm Review

720 s

Frozen Waveform

